

Two PCR Strategy for Locus-Specific Deep Sequencing

Matthew C. Canver¹

¹Harvard Medical School, Boston, MA 02115

1. Extract genomic DNA using Qiagen Blood/Tissue DNeasy kit or other equivalent method.
2. Determine the concentration of each DNA sample using NanoDrop or other equivalent method.
3. Perform PCR #1 reaction(s) using a proofreading enzyme (the following protocol uses Herculase II Fusion DNA Polymerase from Agilent Technologies).
 - Minimize PCR cycles to limit PCR bias. Consider performing multiple different cycle numbers (e.g. 10, 15, 20 cycles) and evaluate on agarose gel.
 - Can try DMSO at increasing concentration 1% to 10%, which often improves specificity. 8% is used in the reaction below.
 - *Amount of genomic DNA (gDNA) for PCR #1 can vary based on experimental needs. On average, a genome from a single cell is approximately 6 picograms. Therefore, 6.6 µg of gDNA represents one million cells. Use adequate gDNA to represent desired number of cells.

PCR #1

X* µL genomic DNA
10 µL of reaction buffer (5x)
1 µL of 100 mM dNTPs
2.5 µL of 5 µM PCR #1 forward primer
2.5 µL of 5 µM PCR #1 reverse primer
4 µL of DMSO
0.5 µL of Herculase II DNA Polymerase
to 50 µL with H₂O

PCR #1 Cycling Conditions

1. 95°C for 2 minutes
2. 95°C for 20 seconds
3. 60°C for 20 seconds
4. 72°C for 30 seconds
5. Repeat steps 2-4 for minimal number of cycles
6. 72°C for 5 minutes

PCR #1 Primers

Forward: **TCGTCCGCAGCGTC**AGATGTGTATAAGAGACAG-Locus-Specific-Sequence

Reverse: **GTCTCGTGGGCTCGG**AGATGTGTATAAGAGACAG-Locus-Specific-Sequence

Blue sequence is Illumina Nextera handle sequence

Recommend 20 bp of locus-specific sequence

4. Perform PCR #2. Each sample will have a unique Illumina Nextera index to allow demultiplexing (see primers below):
 - Minimize cycles to limit PCR bias. Consider performing multiple different cycle numbers (e.g. 10, 15, 20 cycles) and evaluate on agarose gel.

PCR #2

1.0 µL PCR #1 product from step 3 diluted 1:10
2 µL of reaction buffer (5x)
0.1 µL of 100 mM dNTPs
1 µL of 2 µM PCR #2 forward primer
1 µL of 2 µM PCR #2 reverse primer
0.1 µL of Herculase II DNA Polymerase
to 10 µL with H₂O

PCR #2 Cycling Conditions

1. 95°C for 2 minutes
2. 95°C for 20 seconds
3. 60°C for 20 seconds
4. 72°C for 30 seconds
5. Repeat steps 2-4 for minimal number of cycles
6. 72°C for 5 minutes

5. Run the PCR #2 product on an agarose gel and gel purify the band of interest.
6. Quantitate DNA by Qubit or other equivalent method.
7. Perform deep sequencing.

PCR #2 Primers

Forward Primers (**i5-Index-Handle**)

F501 AATGATACGGCGACCACCGAGATCTACAC**TAGATCGCTCGTCGGCAGCGTC**
F502 AATGATACGGCGACCACCGAGATCTACAC**CTCTCTATTCGTCGGCAGCGTC**
F503 AATGATACGGCGACCACCGAGATCTACAC**TATCCTCTTCGTCGGCAGCGTC**
F504 AATGATACGGCGACCACCGAGATCTACAC**AGAGTAGATCGTCGGCAGCGTC**
F505 AATGATACGGCGACCACCGAGATCTACAC**GTAAGGAGTCGTCGGCAGCGTC**
F506 AATGATACGGCGACCACCGAGATCTACAC**ACTGCATATCGTCGGCAGCGTC**
F507 AATGATACGGCGACCACCGAGATCTACAC**AAGGAGTATCGTCGGCAGCGTC**
F508 AATGATACGGCGACCACCGAGATCTACAC**CTAAGCCTTCGTCGGCAGCGTC**
F517 AATGATACGGCGACCACCGAGATCTACAC**GCGTAAGATCGTCGGCAGCGTC**

Reverse Primers (**i7-Index-Handle**)

R701 CAAGCAGAAGACGGCATAACGAGAT**TGCCTTAGTCTCGTGGGCTCGG**
R702 CAAGCAGAAGACGGCATAACGAGAT**CTAGTACGGTCTCGTGGGCTCGG**
R703 CAAGCAGAAGACGGCATAACGAGAT**TTCTGCCTGTCTCGTGGGCTCGG**
R704 CAAGCAGAAGACGGCATAACGAGAT**GCTCAGGAGTCTCGTGGGCTCGG**
R705 CAAGCAGAAGACGGCATAACGAGAT**AGGAGTCCGTCTCGTGGGCTCGG**
R706 CAAGCAGAAGACGGCATAACGAGAT**CATGCCTAGTCTCGTGGGCTCGG**
R707 CAAGCAGAAGACGGCATAACGAGAT**GTAGAGAGTCTCGTGGGCTCGG**
R708 CAAGCAGAAGACGGCATAACGAGAT**CCTCTCTGGTCTCGTGGGCTCGG**
R709 CAAGCAGAAGACGGCATAACGAGAT**AGCGTAGCGTCTCGTGGGCTCGG**
R710 CAAGCAGAAGACGGCATAACGAGAT**CAGCCTCGGTCTCGTGGGCTCGG**
R711 CAAGCAGAAGACGGCATAACGAGAT**TGCCTCTTGTCTCGTGGGCTCGG**
R712 CAAGCAGAAGACGGCATAACGAGAT**TCCTCTACGTCTCGTGGGCTCGG**

Example Locus

Deep sequencing of the -71 DNase hypersensitive site in the *HBS1L-MYB* interval³:

sgRNA (20 bp, chr6:135431630-135431649, hg19)

ACTACTGACATTTATCAACA

PCR #1 primers

Forward: TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCTGCTGGCTTCTTTGCTGTA

Reverse: GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGAGCCTGGGTGACAGAGTGAG

Genomic locus (240 bp, chr6:135431513-135431752, hg19)

CTGCTGGCTTCTTTGCTGTAATATCCTGATCACGCTGACTTCCTTCTGCAACTTCTAGATAAGTAAATTTTTTTT
GATTTATCAGGAAGTGTCTTTGGTCTCTCAGTCAATTCGATTCTACTACTGACATTTATCAACAATGGTGGGTG
TGATATCTTTAATCTAATGAGCTATATAACTGCATTTTCTTTTTGTTGTTGTTGTTGTTTTTTGAGATGGAG
TCTCACTCTGTCACCCAGGCT

PCR #1 amplicon (307 bp amplicon)

TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCTGCTGGCTTCTTTGCTGTAATATCCTGATCACGCTGAC
TTCCTTCTGCAACTTCTAGATAAGTAAATTTTTTTGATTTATCAGGAAGTGTCTTTGGTCTCTCAGTCAATTC
GATTCTACTACTGACATTTATCAACAATGGTGGGTGTGATATCTTTAATCTAATGAGCTATATAACTGCATTT
TCTTTTTGTTGTTGTTGTTGTTTTTTGAGATGGAGTCTCACTCTGTCACCCAGGCTCTGTCTCTTATACACAT
CTCCGAGCCCACGAGAC

PCR #2 amplicon using F501/R701 primers (376 bp amplicon)

AATGATACGGCGACCACCGAGATCTACACTAGATCGCTCGTCGGCAGCGTCAGATGTGTATAAGAGACAGC
TGCTGGCTTCTTTGCTGTAATATCCTGATCACGCTGACTTCCTTCTGCAACTTCTAGATAAGTAAATTTTTTTG
ATTTATCAGGAAGTGTCTTTGGTCTCTCAGTCAATTCGATTCTACTACTGACATTTATCAACAATGGTGGGTGT
GATATCTTTAATCTAATGAGCTATATAACTGCATTTTCTTTTTGTTGTTGTTGTTGTTTTTTGAGATGGAGT
CTCACTCTGTCACCCAGGCTCTGTCTCTTATACACATCTCCGAGCCCACGAGACTAAGGCCGAATCTCGTATGC
CGTCTTCTGCTTG

Blue text = Illumina Nextera handle sequence

Red text = 20 bp of locus specific sequence

Bold text = Illumina Nextera index

Green text = Illumina Nextera adapter

Yellow highlight = sgRNA sequence

Blue highlight = PAM sequence

References

1. Shalem, O. *et al.* Genome-scale CRISPR-Cas9 knockout screening in human cells. *Science* **343**, 84–7 (2014).
2. Canver, M. C. *et al.* BCL11A enhancer dissection by Cas9-mediated in situ saturating mutagenesis. *Nature* **527**, 192–7 (2015).
3. Canver, M. C. *et al.* Variant-aware saturating mutagenesis using multiple nucleases identifies regulatory elements underlying trait-associated DNA. Submitted.
4. Illumina Nextera Deep Sequencing: <http://support.illumina.com>